

PHILIPS

PageWriter

TC20 Cardiograph

**In touch
with you**



Intuitive and efficient

As a world leader in cardiology, we make it our goal to simplify cardiac care, so that you can deliver the best patient care. That is why we created the PageWriter TC20. Advanced, yet easy to use, the PageWriter TC20 offers speed of operation in an attractive and affordable solution that can grow with you as your workflow evolves. Expect fast, efficient clinical workflow, combined with reliable operation for you and your patients.




Advanced, yet easy to use

The PageWriter TC20 cardiograph is designed for use in the demanding hospital environment as well as physician practice settings. The easy 1-2-3 operation, touchscreen display, and additional tools facilitate quality ECG reports.

- Easy 1-2-3 touchscreen operation
- Clinical decision support with the world-class Philips DXL ECG algorithm
- Scalable workflow to meet your needs now and in the future
- Wireless networking capability

The PageWriter TC20 is integrated into a rugged, compact trolley for easy mobility and storage.



Quick and easy placement

Trident 3-in-1 chest leads reduce tangling and lead reversals, for easier placement and quicker ECGs. The anatomical Patient Interface Module mirrors the body, so clinicians can quickly connect the correct leads.

Just touch it

Take ECGs from the touchscreen or the keyboard.

It is as easy as 1-2-3

User-friendly illuminated buttons speed workflow



Connect Leads

The system will perform quality controls, in the form of lead reversal detection and lead checks (impedance).



Enter ID

The ID button enables electronic data entry, reducing the risk of errors caused by entering information manually. Confirm ID with the barcode scanner, or the IECG or EMR interfaces.



Take ECG

Acquire, analyze, print and transmit data with a single button. This standardizes your workflow, so that each ECG is captured and screened, and delivers critical, time-sensitive results to clinicians.



Designed around you

The compact system provides advanced features to support a variety of workflows for you.

Advanced tools at your fingertips


The PageWriter TC20 provides advanced tools to enhance workflow and support clinical decisions. All PageWriter TC cardiographs include the clinical excellence of the DXL ECG algorithm which is built upon over 45 years of research and experience, and provides continuity and consistency in ECG reading and diagnosis through all PageWriter products and the IntelliSpace ECG management system.

Never miss a beat

With Disclosure, review up to five minutes of 12-lead ECG data. Mark events, select pre- and post-event ECGs for analysis, or choose the best traces for use in the ECG report. All at a touch.

Save it

After up to 20 minutes, cardiac event data is automatically saved in a time capsule, so you do not lose sight of an important clinical episode.



IntelliSpace
ECG Management System

Get Orders



PageWriter TC20
Cardiology

Synchronize time

Auto set the PageWriter time with your hospital time master to obtain accurate documentation of your patient's clinical history.

Download orders

Import ECG orders with comprehensive patient information from IntelliSpace ECG (DICOM order manager, EMR).



CCX compliance

Support for open communication protocols like Cisco's CCX and 802.11 (a,b,g,n).


Access ECGs anytime

Ever need a physician to read an ECG 10 minutes after they have left the hospital, or when they are 50 miles away? With TC cardiograph communications capability and IECG, your physicians can access ECGs virtually anywhere for confirmation, over-reads and consultation right around the clock.

New Patient	Worklist 0	Find Patient	Edit Patient ID
Patient ID	123456		
Last Name	Doe		
First Name	Jane		
Gender	Female		
DOB (dd/mm/yyyy)	5/5/1955		
Operator ID	lmg		
Save OK Cancel			

Instant access

Easily acquire or enter patient demographic information by barcode scanning, keyboard entry, worklist download or patient search.




One-touch workflow

The PageWriter TC20 can be configured to automatically print, save and transfer ECG reports – or even retrieve a previous ECG from advanced systems – with a touch of a button.

Clarity when it matters most

Clinical decision support means the right information at the right time presented with clarity to help guide the most productive course of action. The Philips DXL 12-Lead ECG Algorithm provides industry-leading ECG interpretations, particularly with respect to pediatric analysis, pacemaker pulse detection, QT measurements, and a suite of advanced STEMI decision support tools to help guide the treatment of patients with chest pain.




Check and confirm quality

Preview ECG waveforms and interpretation on the 15-inch touchscreen to check for signal quality before printing.

Be sure

Unique LeadCheck software tests for 20 different lead reversals to help you be sure of capturing a diagnostic-quality ECG.



Reveal more

The world-class Philips DXL ECG Algorithm goes beyond traditional 12-lead interpretations to provide STEMI decision support tools and other incremental diagnostic capabilities.

ST map

At a glance, get a clear indication of ST elevation for quick triage.

Critical Values

Quickly identify patients that need urgent care in support of Joint Commission patient safety goals.

STEMI-CA

Culprit artery criteria provide an indication of which artery may be occluded to help you manage your cath lab interventions.

Clinically significant

The previous ECG can be automatically retrieved at the bedside, because a cardiac event is dynamic, with clinical decisions changing frequently during an encounter.

Gender-accurate analysis

Differentiated criteria to help interpret cardiac symptoms in women, including identification of ischemia.

Up-to-date statements

Enhance consistency of care with terminology that conforms to ACC/AHA recommendations.

Higher standards achieved

The solution meets FDA requirements for pharmaceutical studies.

Stay connected

PageWriter TC20 fits seamlessly into your existing IT infrastructure, and supports WiFi 802.11 a/b/g/n. So you are always connected – without being locked in.

Maintain security

The PageWriter TC20 delivers secure, wireless connectivity via standard LAN protocols like 802.11(i) and WPA2/CCMP to protect the privacy of patient, staff, and financial information.



PageWriter TC20 benefits

Cardiology Department Manager

- 1-2-3 buttons light in sequence to guide the user
- User aids support high-quality ECG reports
- High-performance that is cost-effective

Emergency Department Manager

- Quickly and easily see ST elevation
- ECGs needing immediate attention highlighted
- Robust and reliable member of the proven PageWriter TC series

Cardiologist

- Industry-leading Philips DXL ECG Algorithm
- Provides a trusted second opinion
- Advanced STEMI and NSTEMI diagnostic aids

Nurse

- Touchscreen operation is fast and intuitive
- Trident lead wires reduce tangling and lead reversals
- Five minutes of patient ECG history are available for review and creation of ECG reports

IT Administrator

- ECG report export in PDF or XML format
- Advanced wireless capabilities – WPA2, CCX

Technical specifications

ECG functions

Simultaneous lead acquisition	Up to 12 leads
ECG reports: 12-lead	3x4, 3x4 1R, 3x4 3R, 3x4 ST, 3x4 1R ST, 6x2, 12x1, 6x2 1R - Standard and Cabrera formats, plus Pan 12 Cabrera
Rhythm strips	Up to 12 configurable leads
Event marking	Six independent events can be marked for later review and analysis
Full disclosure	- Five minute history of all 12 leads - Complete interpretive ECG report of any 10 seconds
Timed ECG	Support for pharma stress protocols
Report storage/transfer	- Full fidelity at 500 Hz of all 10 sec for all 12 leads - PDF, XML, DICOM 12-lead ECG, and DICOM General ECG formats

Philips DXL 16-lead ECG Algorithm

Interpretive statements	>600 interpretive statements - Integrated pediatric analysis
Leads used in diagnosis	Standard 12 leads plus V3R, V4R, V5R, V7, V8, and V9
Borderline statement suppression	Three configurable settings
Standard measurements	- Ten interval, duration, and axis measurements - Configurable QT correction method
Extended measurements	46 measurements of Morphology analysis in each lead; 21 parameters of Rhythm analysis
Reasons	Selectable explanations of all interpretive statements
Nomenclature	Conforms to 2009 AHA/ACCF/HRS Recommendations for the Standardization and Interpretation of the Electrocardiogram and 2013 ACCF/AHA STEMI Management Guidelines

STEMI diagnostic aids

Graphical ST vector	Two ECG reports with polar ST Maps; frontal and transverse planes
STEMI-CA	Criteria that suggest the probable site of the occlusion
Critical values	Highlights 4 conditions requiring immediate clinical attention

Networked features

Central time management	Time can be manually or automatically synchronized to a network time server via IntelliSpace ECG or IntelliBridge Enterprise
Last ECG orders (requires IntelliSpace ECG)	- Automated retrieval of previous ECG - Configurable rules to retrieve cardiograph-specific Worklists

Signal quality indicators

Leads-off advisory	Anatomical lead map displays the location and label of any loose or disconnected leads/electrodes
Lead color	Four colors to indicate levels of waveform quality
LeadCheck	Lead placement software detects lead reversals
Heart rate	Continuous display of patient heart rate
Print preview	Full screen preview of ECG waveforms to print

Training

Application help	Integrated graphical help screens for primary functions
Self paced	PC-based, interactive, dynamic animation covering all major clinical functionalities

User interface

Touchscreen	1-2-3 operation - Context-sensitive application 5-wire, resistive touchscreen
Keyboard	65-button, standard full alphanumeric keyboard - Special characters supported

Display

Size	6.5 inch TFT active matrix
Resolution	640 x 480 VGA
Colors	64K colors

Signal processing

Sampling rate	8,000 samples per second per lead wire
Patient Interface Module	Remote, microprocessor-controlled digital module provides 5µV resolution

Printer

Resolution	High-resolution, digital-array printer using thermal-sensitive paper; 200 dpi (voltage axis) by 500 dpi (time axis) at 25 mm/s
Paper sizes:	Z-fold letter and A4

Connectivity

LAN connectivity	10/100 Base-T IEEE 802.3 Ethernet via on-board RJ45
Wireless connectivity and security (option D24)	802.11 (a/b/g/n)
Internal storage	200 ECGs
External storage	200 ECGs with optional USB device

Automated data input

Barcode reader (option H17)	- Reads Code 39 Symbolology - Flexible field data entry
Smart "IC" card reader (option H14)	ISO 7816 and EMV 3.1.1; supports SLE 4418/28 and SLE 4443/42

Pre-processing filters

AC noise	50 or 60 Hz
Signal processing	Artifact rejection and baseline wander

Presentation filters – 10 sec reports

High pass	0.05, 0.15, and 0.5 Hz
Low pass	40, 100, and 150 Hz

Presentation filters – rhythm

High pass	0.05 and 0.15 Hz
Low pass	40, 100, and 150 Hz

Electrical

Battery	Lithium ion
Battery capacity	- Typically 30 ECGs on a single charge or 30 minutes of continuous rhythm recording; - No fail operation during ECG printing
Battery recharge	4 hours to full capacity
AC power	100–240 Vac, 50/60 Hz
Power consumption	60 W max

Mechanical

Dimensions	31 x 40 x 21 cm (12 x 16 x 8 in)
Weight	8.6 kg (19 lb) - Includes battery, patient cable

Environmental

Operating conditions	10° to 40°C (50°F to 104°F); 10% to 90% relative humidity (non-condensing); - Up to 4,200 m (14,000 ft) altitude
Storage conditions	(–20°C to 50°C) (–4°F to 122°F); 10% to 90% relative humidity (non-condensing); - Up to 4,550 m (15,000 ft) altitude

Safety and performance

International standards and regulations	- General requirement for safety IEC 60601-1: 1988 +A1:1991 +A2:1995 - Particular requirement for safety of electrocardiographs IEC 60601-2-25: 1993 + A1:1999 - Particular requirements for safety IEC 60601-2-51: 2003 - US general requirements for safety UL 2601-1: 2003 1997 - Diagnostic electrocardiographic devices AAMI EC11 1991 (R: 2001) - CAN/CSA-C22.2 No. 601.1-M90 S1:1994 B:1996 - Electromagnetic compatibility IEC 60601-1-2 second edition 2001
---	--

